

OMNITIP™

Plastic Homogenizer Probes

- **Economical :**
 <1% of the cost of stainless steel probes
- **Eliminates Cross Contamination**
- **7mm (Diameter) x 110mm (Length)**
- **Ideal for PCR, RNA, and DNA sample preparation**
- **Maintenance Free**



Omni Tip™
Plastic Probes

Versatile **Omni Tip™ Probes** are available in both **soft tissue** (deaggregation) and **hard tissue** (frozen) versions. Engineered from durable plastics, Omni Tips™ are reusable (autoclavable) and are so economical they can be disposed of after each use to eliminate cross contamination. The Omni Tips' patented design eliminates the hassle of disassembly and maintenance of traditional stainless steel generators. Treatable volume - 250 µL to 30 mL.

**Gentle
Deaggregation
Omni Tip™**



**Hard Tissue
Omni Tip™**




OMNI
International

TheHomogenizerCompany.com

OMNI TIP™

Plastic Homogenizer Probes

Omni Tip™ Plastic Probes are available in the convenient **Omni Tip™ Homogenizing Kits**. They include 12 Omni Tips™, the reliable Omni TH-Tissue Homogenizer, Omni Tip™ Adapter, and an optional stainless steel probe, all in a convenient carrying case.



Omni Tip™ Plastic Probes are also used with the revolutionary **Omni Prep Multi Sample Homogenizer**. Homogenize six samples with the touch of a button!



Omni Tip™ Customer Testimonials

With Omni Tips I can be certain that all the 'tissue sample has been homogenized; there is no sample loss. I found that the Omni Tip hard tissue version worked as well as the stainless steel probe."

- Rob T.
National Institutes of Health

"[The Omni Tips] worked very well. Will use these in the future. Much better than the method I currently use. Thanks!"

- Thud S.
Cornell University

"The hard tissue Omni Tips produce great yield; 10 times greater yield when compared to a conventional mortar and pestle!"

- Fernando N.
Arcadia Biosciences

"...We use two Omni homogenizers and have found the disposable tip plastic probes work great and we don't get cross contamination between samples. For mouse samples we use the standard tips and for tougher tissues, the hard tissue plastic probe tips. Both work very well."

- Julie E.
UPMC

"[The Hard Tissue Tips] work much better than traditional blades for cutting through muscle, brown adipose tissue and uterus. Blades are also durable enough for several uses."

- Alison B.
University of Massachusetts, Amherst

"[The Hard Tissue Omni Tips] work just as well as stainless steel probes; no worries about maintenance of probe."

- Debbie C.
University of Pennsylvania

The Omni hard tissue tips have been indispensable to us in the analysis of heart, skeletal muscle and tumor tissue. Thank you for this product.

- Caryl A.
Methodist Research Center



935-C Cobb Place Blvd
Kennesaw, GA 30144
omni@omni-inc.com

800-776-4431
770-421-0058
770-421-0206 Fax

www.omni-inc.com