



Using CTSV FILCONS with different mesh sizes, you can easily remove cell clumps and large particles from your sample after tissue dissociation and ensure that you have a uniform single-cell suspension.

The strainers have smart features, such as:

- Compatibility with 5 (FILCONS) and 50 ml (FILCONS MAX) conical tubes
- Prevent clogging with optimized ventilation
- Compact size for sterile handling
- Syringe fitting for bigger volumes, cup type for smaller volumes

FILCONS 20 μ m



020-22N-S



120-22N-S



WWW.CTSV.BIZ



FILCONS 20 μ m

APPLICATION TABLE

20 μ m

STRAINER TARGET

Leydig and Sertoli cells, sperm, dendritic cells, monocytes

TISSUES


Testicle, Lymph Nodes, Blood, Intestine

DISSOCIATION SYSTEM/DEVICE

Medicons-P cod. 79700 (70 μ m cut)



Testicle



WWW.CTSV.BIZ