



FILCONS 200µm

Using CTSV FILCONS with different mesh sizes, you can easily remove cell clumps and large particles from your sample after tissue dissociation and ensure that you have a uniform single-cell suspension.

The strainers have smart features, such as:

- Compatibility with 5 (FILCONS) and 50 ml (FILCONS MAX) conical tubes
- Prevent clogging with optimized ventilation
- Compact size for sterile handling
- Syringe fitting for bigger volumes, cup type for smaller volumes



200-20N-S

220-20N-S



WWW.CTSV.BIZ

FILCONS 200µm

APPLICATION TABLE

	200 µm
STRAINER TARGET	Oocytes, parenchyma cells, Protozoa, some species of Euglena
TISSUES	Ovary, Various plant tissues (parenchyma, leaves, stems, roots), Aquatic environments
DISSOCIATION SYSTEM/DEVICE	Medicons cod.79400 (50 µm cut) Medicons Max cod. 79500 for plant tissues



WWW.CTSV.BIZ