



FILCONS 500 μ m

Using CTSV FILCONS with different mesh sizes, you can easily remove cell clumps and large particles from your sample after tissue dissociation and ensure that you have a uniform single-cell suspension.

The strainers have smart features, such as:

- Compatibility with 5 (FILCONS) and 50 ml (FILCONS MAX) conical tubes
- Prevent clogging with optimized ventilation
- Compact size for sterile handling
- Syringe fitting for bigger volumes, cup type for smaller volumes



500-50N-S

550-50N-S



WWW.CTSV.BIZ



APPLICATION TABLE

	500 μm
STRAINER TARGET	Macroscopic algae such as brown algae or red algae, some insect eggs, some mushroom spores, some plant ovules such as conifer ovules
TISSUES	Land or aquatic environment
DISSOCIATION SYSTEM/DEVICE	-



WWW.CTSV.BIZ